

Activity 5

Aim: To detect the presence of adulterant in given food samples.

Theory: The substance which lowers or degrades the quality of food material is called adulterant.

Sr. No.	Name of food material	Common Adulterant
1	Powdered Sugar	Baking soda
2	Turmeric powder	Metanil yellow color
3	Tea powder	Exhausted tea leaves dried and artificially Colored
4	Chili powder	Brick powder
5	Edible oil	Argemone oil

Chemicals: Given food samples, dil. HCl, conc. HCl, etc.

Apparatus: Test tubes, test tube stand, etc.

Procedure:

Sr. No.	Name of Food	Test	Observation	Inference
1	Powdered sugar	Sample of sugar + dil. HCl	Effervescence	Sugar sample is adulterated
			No effervescence	Sugar sample is not adulterated / Pure
2	Turmeric powder	Sample of turmeric powder + few drops of conc. HCl	Violet color	Turmeric powder sample is adulterated
			No violet color	Turmeric powder sample is not adulterated / Pure
3	Tea powder	Sprinkle the sample of tea powder on a wet filter paper	Yellow or red colored spot appears	Tea sample is adulterated
			No yellow or red spot appears	Tea sample is not adulterated / Pure
4	Chili powder	Sprinkle chili powder on the surface of water	Brick powder settled at the bottom of test tube	Chili powder sample is adulterated
			No brick powder settled at the bottom of test tube	Chili powder sample is not adulterated / Pure
5	Edible oil	5 ml sample of edible oil + 5 ml of conc. HCl and shake.	Yellow or orange yellow color appears in lower acid layer	Edible oil sample is adulterated
			No yellow or orange yellow color appears in lower acid layer	Edible oil sample is not adulterated / Pure

